

## LZCap®AG(3'Ma-Cy7) Capping Kit (N1-Me-pUTP)

**Description:** LZCap®AG(3'Ma-Cy7) Capping Kit is the co-transcriptional synthesis of mRNA containing the Cap1 structure *in vitro*. The core component LZCap®AG(3'Ma-Cy7) is a Cap1 analog that can be added to the 5' end of mRNA in the presence of T7 polymerase. DNase I and LiCl are included in the kit for mRNA purification. Recommended wavelengths for detecting the Cy7 fluorescence of LZCap®AG(3'Ma-Cy7) mRNA: (651/780).


**Specifications and Components:** 20T/Kit (20µL)

No.	Component	Concentration	20 Test
1	LZCap®AG(3'Ma-Cy7)	25 mM	64 µL
2	ATP	100 mM	20 µL
3	N1-Me-pUTP	100 mM	20 µL
4	CTP	100 mM	20 µL
5	GTP	100 mM	20 µL
6	10×Fluro Transcription Buffer	/	40 µL
7	Enzyme Mix	/	80 µL
8	Recombinant DNase I(RNase-free)	5 U/µL	40 µL
9	LiCl	7.5 M	1 mL
10	Control	0.5 µg/µL	4 µL
11	RNase Free Water	/	2 mL

**Storage Conditions:** store at -15°C or below.

### DNA Template Design

LZCap®AG(3'Ma-Cy7) is suitable for AG-initiated sequences. As shown in the figure below, the T7 promoter (underlined) followed by the AG sequence can effectively initiate transcription.



## Protocol

### 1. Thaw components required for the experiment on ice.

Component	Volume (μL)	Final concentration
ATP (100 mM)	1μL	5 mM
N1-Me-pUTP (100 mM)	1μL	5 mM
CTP (100 mM)	1μL	5 mM
GTP (100 mM)	1μL	5 mM
LZCap®AG(3'Ma-Cy7) (25 mM)	3.2μL	4 mM
10×Fluro Transcription Buffer	2μL	1×
Enzyme Mix	4μL	/
Linear DNA+ RNase Free Water	6.8μL	50 ng/μL
终体积	20μL	

- Mix the prepared reaction solution, centrifuge briefly, and incubate at 37°C for 2-3 hours. The reaction time should be increased to 4-8 h when the length of transcript is less than 100 nt.
- After the reaction, 2μL of DNase I is added to each tube and the DNA template will be digested at 37°C for 15min.
- Purification of mRNA by LiCl precipitation
  - Add 50μL of LiCl and 30μL of RNase Free Water to 10μL of transcript mRNA after the reaction (the final concentration of LiCl should be kept at 4.5-4.8 M), mix well and incubate at -20°C for at least 0.5h.
  - Centrifuge the mixture at 12000rpm for 15min, and remove the supernatant and preserve the precipitate.
  - Wash the precipitate with 600μL of pre-chilled 75% ethanol, centrifuge at 12000 rpm for 8 min, and remove the supernatant.
  - The purified mRNA should be dried for 10 min until the ethanol evaporated completely and is re-dissolved with 30-100μL RNase-Free water.

## Notes:

- LZCap®AG(3'Ma-Cy7) is suitable for T7 promoter transcription vector with 5' AG 3' initiated sequences, which needs to be considered when constructing the vector.
- LZCap®AG(3'Ma-Cy7) and its mRNA products should be stored and used away from light.
- The reagents, consumables and containers used in this experiment are free of RNase and DNase contamination.
- It is recommended to use a linearized DNA template for transcription.

- 5) When modified nucleotides is used in place of wild-type nucleotides, the final concentration of the transcript is not affected, but the UV ratio of 260/280 may be different from regular RNA.
- 6) If the PCR product is used as the DNA template, the amount of DNA template can be reduced by half.
- 7) Due to the high concentration of 10×Fluro Transcription Buffer, the high salt environment will lead to polymerase inactivation. When preparing the reaction solution, we need to add water first, then buffer, NTPs and LZCap®AG(3'Ma-Cy7), DNA template, and finally the enzymes.